



MD2030B

REF	PART NO.	QTY	DESCRIPTION	NOTES
1	3580382	1	Injection pump	
2		1	• Housing	
3	3580402	3	• Delivery valve retai	
4	3580124	3	• O-ring	
5	3580125	3	• Spring	
6	3580276	3	• Gasket	
7	(3580128)	3	• Delivery valve	Obsolete part
8	3580403	3	• Delivery valve	
9	3580130	3	• Control sleeve	
10	3580131	3	• Washer	
11	3580132	3	• Spring	
12	3580133	3	• Washer	
13	(3580134)	X	• Adjusting washer	TH=1,85mm, Obsolete part
	(3580135)	X	• Adjusting washer	Obsolete part, TH=1,90mm
	(3580136)	X	• Adjusting washer	Obsolete part, TH=1,95mm
	(3580137)	X	• Adjusting washer	Obsolete part, TH=2.00mm
	(3580138)	X	• Adjusting washer	Obsolete part, TH=2.05mm
	(3580139)	X	• Adjusting washer	Obsolete part, TH=2.10mm
	(3580140)	X	• Adjusting washer	Obsolete part, TH=2.15mm
	(3580141)	X	• Adjusting washer	Obsolete part, TH=2.20mm
	(3580142)	X	• Adjusting washer	Obsolete part, TH=2.25mm
	(3580143)	X	• Adjusting washer	Obsolete part, TH=2.30mm
	(3580144)	X	• Adjusting washer	Obsolete part, TH=2.35mm
	3580145	X	• Adjusting washer	TH=2.40mm
	3580146	X	• Adjusting washer	TH=2.45mm
	3580147	X	• Adjusting washer	TH=2.50mm
14	3580148	3	• Tappet	
15	3580278	1	• Control rod	
16		1	• Bracket	
17	3580151	2	• Screw	
18	3580152	1	• Screw	
19	18665	2	• Gasket	
20	3580154	1	• Banjo nipple	
21	3580155	3	• Lock plate	
22	3580156	3	• Screw	
23	3580157	3	• Locking pin	